REF 13-09A-100

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- Methods for measurement of gastric motility
Lawrence A. E., *et al.*
Am J Physiol Gastrointest Liver Physiol (2009) 296, 461-75
- Measurement of gastric emptying rate of solids by means of a carbon labeled octanoic acid breath test
Ghoos Y.F., *et al.*
Gastroenterology (1993) 104, 1640-47
- Gastric emptying and intragastric balloon in obese patients
Bonazzi P., *et al.*
European Review for Medical and Pharmacological Sciences (2005) 9(Suppl 1), 15-21
- 13C octanoic acid breath test for measuring gastric emptying of solids
Perri F., *et al.*
Eur. Rev. Med. Pharmacol. Sci. (2005) 9, 3-8
- The Wagner-Nelson method makes the [13C]-Breath Test comparable to radioscintigraphy in measuring gastric emptying of solid / liquid test meal in human
Sanaka M. *et al.*
Clinical and Experimental Pharmacology and Physiology (2007) 34, 641-44
- 13C-Octanoic acid breath test for gastric emptying rate of solids
Maes B., *et al.*
Gastroenterology (1998) 114, 856-63
- Procedures in the 13C octanoic acid breath test for measurement of gastric emptying: analysis using Bland-Altman methods
Clegg M.E., *et al.*
Scandinavian Journal of Gastroenterology (2010) 45, 852-61
- Gastric emptying flow curves separated from carbon-labeled octanoic acid breath test results
Maes B., *et al.*
Am J Physiol Gastrointest Liver Physiol (1998) 275, 169-75
- Scintigraphic evaluation of gastric emptying
Rehm P.
Applied Radiology (2002) 26-32
- Long term effects of lipase inhibition by orlistat on gastric emptying and OCT for a solid meal
Kocelak P., *et al.*
J Gastroenterol (2008) 43, 609-17
- Delayed gastric emptying in the obese: an assessment using thenon-invasive 13C-octanoic acid breath test
Jackson S. J., *et al.*
Diabetes, Obesity and Metabolism (2004) 6, 264-70
- 13C breath test for patients with non alcoholic steatohepatitis
Mawatari H., *et al.*
Hepato Gastroenterologi (2009) 56, 1436-38
- Pharmacological modulation of gastric emptying rate of solids as measured by the carbon labeled octanoic acid breath test.
Maes B., *et al.*
Gut (1994) 35, 333-37
- Influence of cisapride on gastric emptying of solids and liquids monitored by the 13C breath test
Duan L., *et al.*
Dig. Dis. Sci. (1995), 50, 2200-06
- Influence of octreotide on gastric emptying of solids and liquids in normal healthy subjects.
Maes B., *et al.*
Aliment. Pharmacol. Ther. (1995) 9, 11-18
- Application of isotope-selective nondispersive infrared spectrometry (IRIS) for evaluation of [13C]octanoic acid gastric-emptying breath tests: comparison with isotope ratio-mass spectrometry (IRMS)
Schadewaldt P., *et al.*
Clinical Chemistry (1997) 43:3, 518-22
- The accuracy of NDIRS vs IRMS
Leodolter A.
Gastroenterology 118 no 4
- 13C octanoic acid breath test for non invasive assessment of gastric emptying in diabetic patients
Ziegler D., *et al.*
Diabetologia (1996) 39, 823-30
- 13C octanoic acid breath test for gastric emptying measurements
Delbende B., *et al.*
Eur. J. Gastroenterol Hepatol (2000) 12, 85-91
- Early satiety is the only patient reported symptom associated with gastric emptying, as assessed by breath test
Ron Y., *et al.*
J. Neuro Gastroenterol Motil (2011) 17, 61-66
- A valid, accurate, office based non-radioactive test for gastric emptying of solids
Lee J.S., *et al.*
Gut (2000) 46, 768-73
- Gastric emptying of fluids and solids in healthy subjects
Hellmig S., *et al.*
J GastroenterolHepatol (2006) 12, 1832-8
- Gastric emptying in diabetic patients by the 13C-octanoic acid breath test: role of insulin in gastric motility
M. Matsumoto *et al.*
J Gastroenterol (2007); 42:469-474
- Two point, time saving method for measurement of gastric emptying with diagnostic accuracy compared to that of the conventional method
Halkar R.
Radiology (1999) 213, 599-602

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